



DISPOSABLE PINS FOR CRANIAL FIXATION



FEATURES



Tip of the adult model pin.



Tip of the pediatric model pin.

- Angled tip with great mechanical resistance. Its great angle prevents torsion or the tip bending against the patient's bone when applying force;
- Its design that offers great fixation and minimum impact on the patient;

- Ergonomic that facilitates the insertion and removal of the pin in the cranial fixator;
- Adult and pediatric models, especially designed to be minimally invasive;
- Sterilization by ETO.


SECURITY

The Micromar cranial fixation pins have a protection at the tip to avoid accidents during the patient preparation.



MODELS

Micromar offers pins for cranial fixation produced in Stainless Steel, Titanium and Aluminum - the latter lighter and with less artifact interference in the images - developed with different characteristics to meet different surgical needs.



The single-use cranial fixation pin kit consists of 03 units of the same model.



CODE	MATERIAL	USE
MD0300FAA	Stainless Steel	Pediatric
MD0100FAA	Stainless Steel	Adult



MD0310FAA	Aluminum	Pediatric
MD0110FAA	Aluminum	Adult



MD0320FAA	Titanium	Pediatric
MD0120FAA	Titanium	Adult

MICROMAR

INSPIRED BY LIFE



ING - VER: 03 - 01/06/2023. All rights reserved.

All images shown in this catalog are for illustration purposes only. Check the current version of the brochure available on our website.

Micromar Ind. Com. Ltda.

Av. Marginal ao Córrego da Serraria, nº 168

CEP: 09980-390


Diadema · São Paulo · Brasil

+55 11 4057-5722

international@micromar.com

 micromar.com

 facebook.com/micromar

 instagram.com/micromarneurosurgery

 linkedin.com/company/micromar