


# EasyDrill®

Cranial Perforator



M I C R O M A R  
INSPIRED BY LIFE



**CE | FDA\***

**4 YEARS VALIDITY**

## EasyDrill® Efficient cutting combined with high technology.

Based on an intense research, Micromar has used all its tradition in neurosurgery to develop the cranial perforator for burr hole trepanation EasyDrill®.

Designed to be used with Hudson couplings, provides benefits such as:

- Autostop technology;
- Exclusive cut geometry – quick and smooth perforations;
- Non-slip tip – stability and accuracy during drilling;
- Pre-tested one by one;
- Worldwide distributed by Medtronic®.



Safe non-slip start prevents EasyDrill® from sliding on the skull.

Stable from the start!



### Dura protection

Its unique, patented cammed-lug release mechanism, has been developed to meet the strictest quality standards and rigorous neurosurgeons needs.

EasyDrill's auto-stop mechanism minimizes the risk of Dura's dangerous penetration, allowing a faster and more reliable drilling.



Code	BA0008FAA	BA0011FAA	BA0210FAA	BA0010FAA
Outer/Inner Diameter	9/6 Ø	11/7 Ø	14/11 Ø	14/11 Ø
Bone Shelf	 3 mm	 3 mm	 1,5 mm	 3 mm

# Quality and storage control



EasyDrill® Cranial Perforator Micromar have a strict quality control, in which each drill is pre-tested before being packaged and made available for surgery.

- Dedicated package;
- Complete manual following International regulations;
- Already sterilized – saving time in the surgical procedure;
- Labeling with barcode to facilitate inventory control;
- Total traceability of the products in the field.

Main Distributor:






M I C R O M A R  
INSPIRED BY LIFE

ING - VER: 08 - 08/05/2024. All rights reserved.  
All images shown in this catalog are for illustration purposes only. Current version of the brochure available on our website.  
Check availability and the official distributor in your region. ANVISA REGISTER: 80051250022

**Micromar Ind. Com. Ltda.**

Av. Marginal ao Córrego da Serraria, nº 168  
CEP: 09980-390  
Diadema · São Paulo · Brasil  
+55 11 4057-5722  
international@micromar.com

 micromar.com  
 facebook.com/micromar  
 instagram.com/micromarneurosurgery  
 linkedin.com/company/micromar